



Anti-LASV GPC Monoclonal Antibody, 25.6A

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT REFERENCE | PTX19115 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1, kappa |
| TARGET | Pre-glycoprotein polyprotein GP complex, Pre-GP-C, GPC, GP-C, Lassa Virus (LASV) |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Human |
| REACTIVITY | Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV) |
| FORM & BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

