



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-LASV GPC Monoclonal Antibody, 36.9F

OVERVIEW

PRODUCT REFERENCE	PTX19117
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Pre-glycoprotein polyprotein GP complex, Pre-GP-C, GPC, GP-C, Lassa Virus (LASV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.