



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-HRSV-A G/Major surface glycoprotein G Monoclonal Antibody, 2D10

OVERVIEW

PRODUCT REFERENCE	PTX19122
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	G, Major surface glycoprotein G, Attachment glycoprotein G, Membrane-bound glycoprotein, mG, Human Respiratory Syncytial Virus (HRSV), RSV

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human respiratory syncytial virus A (strain A2)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.