



Anti-HRSV-A G/Major surface glycoprotein G Monoclonal Antibody, 3D3

OVERVIEW

PRODUCT REFERENCE	PTX19123
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	G, Major surface glycoprotein G, Attachment glycoprotein G, Membrane-bound glycoprotein, mG, Human Respiratory Syncytial Virus (HRSV), RSV

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human respiratory syncytial virus A (strain A2)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

