



## Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Monoclonal Antibody, AR3C

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19130
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein, Hepatitis C Virus (HCV)

### PRODUCT INFORMATION

---

<b>HOST</b>	Human
<b>REACTIVITY</b>	Hepatitis C virus genotype 1a (isolate H77) (HCV)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

