



Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Monoclonal Antibody, HC84.26.5D

OVERVIEW

PRODUCT REFERENCE	PTX19131
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein, Hepatitis C Virus (HCV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Hepatitis C virus genotype 1a (isolate H77) (HCV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

