



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Monoclonal Antibody, C2

OVERVIEW

PRODUCT REFERENCE	PTX19133
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein, Hepatitis C Virus (HCV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Hepatitis C virus genotype 1a (isolate H77) (HCV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.