



## Anti-CHIKV Spike glycoprotein E2 Monoclonal Antibody, 5M16

### OVERVIEW

---

|                          |  |
|--------------------------|--|
| <b>PRODUCT REFERENCE</b> | PTX19138   |
| <b>CLONALITY</b>         | Monoclonal Antibody  |
| <b>ISOTYPE</b>           | IgG1, kappa  |
| <b>TARGET</b>            | Structural polyprotein, Envelope glycoprotein E2, Spike glycoprotein E2, Chikungunya Virus (CHIKV) |

### PRODUCT INFORMATION

---

|                          |   |
|--------------------------|---|
| <b>HOST</b>              | Human   |
| <b>REACTIVITY</b>        | Chikungunya virus (strain S27-African prototype) (CHIKV)  |
| <b>FORM &amp; BUFFER</b> | 0.01M PBS, pH 7.4.  |
| <b>STORAGE</b>           | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

