



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-CHIKV Spike glycoprotein E2 Monoclonal Antibody, 5F10

OVERVIEW

PRODUCT REFERENCE	PTX19140
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Structural polyprotein, Envelope glycoprotein E2, Spike glycoprotein E2, Chikungunya Virus (CHIKV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Chikungunya virus (strain S27-African prototype) (CHIKV), Chikungunya virus (strain Nagpur) (CHIKV), Chikungunya virus (strain 37997) (CHIKV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.