



Anti-Mouse CD274/PD-L1/B7-H1 Monoclonal Antibody, SAA0142

OVERVIEW

PRODUCT REFERENCE	PTX19147
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-Kappa
TARGET	B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

