



## Anti-Human ADAMTS5 Monoclonal Antibody, 12F4

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19150
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	ADAMTS5, ADAM-TS 5, ADAM-TS5, ADAMTS-5, ADAMTS11, ADAM-TS 11, ADAMTS-11, ADMP2, ADMP-2, A disintegrin and metalloproteinase with thrombospondin motifs 5, A disintegrin and metalloproteinase with thrombospondin motifs 11, Aggrecanase-2

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

