



Anti-Human ADAMTS5 Monoclonal Antibody, 7B4

OVERVIEW

PRODUCT REFERENCE	PTX19151
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	ADAMTS5, ADAM-TS 5, ADAM-TS5, ADAMTS-5, ADAMTS11, ADAM-TS 11, ADAMTS-11, ADMP2, ADMP-2, A disintegrin and metalloproteinase with thrombospondin motifs 5, A disintegrin and metalloproteinase with thrombospondin motifs 11, Aggrecanase-2

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

