



Anti-Human ADAMTS4 Monoclonal Antibody, 7C7

OVERVIEW

PRODUCT REFERENCE	PTX19153
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	ADAMTS4, A disintegrin and metalloproteinase with thrombospondin motifs 4, ADAM-TS 4, ADAM-TS4, ADAMTS-4, ADMP-1, Aggrecanase-1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

