



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human FALZ/BPTF IP Monoclonal Antibody, SAA0188

OVERVIEW

PRODUCT REFERENCE	PTX19156
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	Fetal Alzheimer antigen, FALZ, Nucleosome-remodeling factor subunit BPTF, FAC1, BPTF, Fetal Alz-50 clone 1 protein, Bromodomain and PHD finger-containing transcription factor

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human, Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.