



Anti-Human CIAO1 IP Monoclonal Antibody, SAA0198

OVERVIEW

PRODUCT REFERENCE	PTX19157
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	WD repeat-containing protein 39, Probable cytosolic iron-sulfur protein assembly protein CIAO1, WDR39, CIA1, CIAO1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

