



## Anti-Human USP7 IP Monoclonal Antibody, SAA0238

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19164
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	Ubiquitin-specific-processing protease 7, Deubiquitinating enzyme 7, Ubiquitin thioesterase 7, USP7, Herpesvirus-associated ubiquitin-specific protease, HAUSP, Ubiquitin carboxyl-terminal hydrolase 7

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

