



Anti-Human ASF1B IP Monoclonal Antibody, SAA0242

OVERVIEW

PRODUCT REFERENCE	PTX19165
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	ASF1B, hAsf1, hCIA-II, CIA-II, hAsf1b, Histone chaperone ASF1B, Anti-silencing function protein 1 homolog B, CCG1-interacting factor A-II

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

