



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human TRIM28 IP Monoclonal Antibody, SAA0249

OVERVIEW

PRODUCT REFERENCE	PTX19167
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	TIF1-beta, Tripartite motif-containing protein 28, KRAB-associated protein 1, TRIM28, RING-type E3 ubiquitin transferase TIF1-beta, KRAB-interacting protein 1, KRIP-1, E3 SUMO-protein ligase TRIM28, RNF96, Nuclear corepressor KAP-1, TIF1B, Transcription i

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.