



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human TIF1/TRIM24 IP Monoclonal Antibody, SAA0265

OVERVIEW

PRODUCT REFERENCE	PTX19168
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	RING-type E3 ubiquitin transferase TIF1-alpha, TIF1-alpha, RING finger protein 82, Transcription intermediary factor 1-alpha, RNF82, Tripartite motif-containing protein 24, E3 ubiquitin-protein ligase TRIM24, TIF1, TIF1A, TRIM24

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.