



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human TROP2/TACSTD2/TROP-2 Monoclonal Antibody, FU-MK1

OVERVIEW

PRODUCT REFERENCE	PTX19170
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Pancreatic carcinoma marker protein GA733-1, Tumor-associated calcium signal transducer 2, M1S1, TROP2, TACSTD2, Cell surface glycoprotein Trop-2, GA733-1, Membrane component chromosome 1 surface marker 1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.