



Anti-Mouse PCSK9 Monoclonal Antibody, mAb1

OVERVIEW

PRODUCT REFERENCE	PTX19171
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Proprotein convertase subtilisin/kexin type 9, Subtilisin/kexin-like protease PC9, NARC-1, PC9, Proprotein convertase 9, Neural apoptosis-regulated convertase 1, NARC1, PCSK9

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Human, Mouse, Monkey
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

