



Anti-Mouse MMP9 Monoclonal Antibody, L51/82

OVERVIEW

PRODUCT REFERENCE	PTX19175
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Matrix metalloproteinase-9, CLG4B, GELB, 92 kDa type IV collagenase, 92 kDa gelatinase, Gelatinase B, MMP9, MMP-9

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human, Mouse, Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

