



Anti-Mouse GPC3 Monoclonal Antibody, 7C8

OVERVIEW

PRODUCT REFERENCE	PTX19179
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	GPC3, MXR7, GTR2-2, Glypican-3, Intestinal protein OCI-5, OCI5

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

