



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Mouse SERPINE1/PAI-1 Monoclonal Antibody, SAA0768

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT REFERENCE | PTX19180 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1, kappa |
| TARGET | PLANH1, Endothelial plasminogen activator inhibitor, SERPINE1, Plasminogen activator inhibitor 1, PAI1, PAI, PAI-1, Serpin E1 |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Human |
| REACTIVITY | Mouse, Rat |
| FORM & BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.