



Anti-Mouse STAT3 Monoclonal Antibody, SAA0274

OVERVIEW

PRODUCT REFERENCE	PTX19183
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Signal transducer and activator of transcription 3, APRF, STAT3, Acute-phase response factor

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Human, Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

