



Anti-Mouse TRPC1 Monoclonal Antibody, 1F1R

OVERVIEW

PRODUCT REFERENCE	PTX19185
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	TRPC1, TRP-1, TRP1, TrpC1, Short transient receptor potential channel 1, Transient receptor protein 1

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Human, Mouse, Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

