



Anti-Mouse IgE/IGHE Monoclonal Antibody, 84.1C

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT REFERENCE | PTX19186 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1, kappa |
| TARGET | Ig epsilon chain C region, IGHE, Ig epsilon chain C region ND, Immunoglobulin heavy constant epsilon |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Rat |
| REACTIVITY | Mouse |
| FORM & BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

