



## Anti-Mouse CXCL12/SDF-1 Monoclonal Antibody, 30D8

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19187
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	Intercrine reduced in hepatomas, hSDF-1, Stromal cell-derived factor 1, SDF-1, Pre-B cell growth-stimulating factor, SDF1, SDF1A, C-X-C motif chemokine 12, SDF1B, CXCL12, IRH, PBSF, hIRH

### PRODUCT INFORMATION

---

<b>HOST</b>	Rat
<b>REACTIVITY</b>	Human, Mouse
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

