



Anti-Mouse CXCL10/IP-10 Monoclonal Antibody, 54

OVERVIEW

PRODUCT REFERENCE	PTX19188
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Gamma-IP10, C-X-C motif chemokine 10, INP10, 10 kDa interferon gamma-induced protein, Small-inducible cytokine B10, SCYB10, IP-10, CXCL10

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Human, Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

