



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Mouse IL2 Monoclonal Antibody, JES6-1A12

OVERVIEW

PRODUCT REFERENCE	PTX19190
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b, kappa
TARGET	IL2, IL-2, Aldesleukin, Interleukin-2, TCGF, T-cell growth factor

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.