



Anti-Mouse IL4 Monoclonal Antibody, 11B11

OVERVIEW

PRODUCT REFERENCE	PTX19192
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Binetrakin, B-cell stimulatory factor 1, Pitrakinra, Lymphocyte stimulatory factor 1, IL-4, BSF-1, IL4, Interleukin-4

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

