



Anti-Mouse IL22 Monoclonal Antibody, 368D04

OVERVIEW

PRODUCT REFERENCE	PTX19195
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Interleukin-22, Cytokine Zcyto18, IL-22, ZCYTO18, IL22, ILTIF, IL-TIF, IL-10-related T-cell-derived-inducible factor

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Human, Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

