



Anti-Human IL25/IL17E Monoclonal Antibody, 2c3

OVERVIEW

PRODUCT REFERENCE	PTX19196
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Interleukin-25, IL-25, IL17E, IL25, IL-17E, Interleukin-17E

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

