



Anti-Human TYRP1 Monoclonal Antibody, TA99

OVERVIEW

PRODUCT REFERENCE	PTX19205
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	5, 6-dihydroxyindole-2-carboxylic acid oxidase, Glycoprotein 75, CAS2, TYRP1, TRP-1, TYRRP, TRP1, Catalase B, TYRP, Melanoma antigen gp75, TRP, DHICA oxidase, Tyrosinase-related protein 1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human, Mouse, Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

