



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Phospho-AKT1 (pS212) Monoclonal Antibody, SAA0279

OVERVIEW

PRODUCT REFERENCE	PTX19208
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	RAC-alpha serine/threonine-protein kinase, RAC, RAC-PK-alpha, PKB alpha, Proto-oncogene c-Akt, Protein kinase B alpha, Protein kinase B, AKT1, PKB

PRODUCT INFORMATION

HOST	Rat
REACTIVITY	Human, Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.