



Anti HLA-G Monoclonal Antibody, 15E7#

OVERVIEW

PRODUCT REFERENCE	PTX19209
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	HLA-G, HLAG, sHLA-G, HLA G antigen, MHC class I antigen G, HLA-6.0, HLA class I histocompatibility antigen, alpha chain G, HLA-G-0031

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

