



Anti-Human Transthyretin amyloidosis Monoclonal Antibody, 5A1

OVERVIEW

PRODUCT REFERENCE	PTX19213
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b, kappa
TARGET	Transthyretin, Prealbumin, TBPA, TTR, PALB, ATTR, Transthyretin amyloidosis, mis-TTR

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

