



Anti-Vibrio cholerae ctxB/Cholera Toxin Subunit B Monoclonal Antibody, TE33

OVERVIEW

PRODUCT REFERENCE	PTX19215
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Cholera enterotoxin subunit B, Cholera enterotoxin B chain, Cholera enterotoxin gamma chain, Cholera toxin B subunit, CTxB, CTB, VC_1456, Bacteria

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Vibrio cholerae serotype O1 (strain ATCC 39315 / El Tor Inaba N16961), Escherichia coli
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

