



Anti-Mouse KIRREL3/NEPH2 Monoclonal Antibody, N320/48

OVERVIEW

PRODUCT REFERENCE	PTX19217
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	KIRREL3, NEPH2, Kin of irregular chiasm-like protein 3, Nephrin-like protein 2, KIAA1867, Kin of IRRE-like protein 3

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

