



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-BVDV Envelope glycoprotein E2/gp55 Monoclonal Antibody, H15

OVERVIEW

PRODUCT REFERENCE	PTX19218
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Genome polyprotein, Envelope glycoprotein E2, gp55, Bovine viral diarrhea virus (BVDV)

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Bovine viral diarrhea virus (BVDV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.