



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Phospho-Smad3 (pS208/pS213) Monoclonal Antibody, 5A11

OVERVIEW

PRODUCT REFERENCE	PTX19220
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Mad3, JV15-2, Mothers against decapentaplegic homolog 3, SMAD family member 3, Smad3, SMAD3, MADH3, hMAD-3, SMAD 3, MAD homolog 3, hSMAD3, Mothers against DPP homolog 3

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.