



Anti-Phospho-RIP3 (pS227) Monoclonal Antibody, P2.A11

OVERVIEW

PRODUCT REFERENCE	PTX19231
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	RIPK3, RIP-like protein kinase 3, RIP3, Receptor-interacting serine/threonine-protein kinase 3, RIP-3, Receptor-interacting protein 3

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

