



Anti-Phospho-RIP3 (pS202) Monoclonal Antibody, P3.28

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT REFERENCE | PTX19232 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG2a, kappa |
| TARGET | RIPK3, RIP-like protein kinase 3, RIP3, Receptor-interacting serine/threonine-protein kinase 3, RIP-3, Receptor-interacting protein 3 |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Mouse |
| REACTIVITY | Human |
| FORM & BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

