



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Phospho-HTRA2 (pS212) Monoclonal Antibody, P6.C12

OVERVIEW

PRODUCT REFERENCE	PTX19237
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	High temperature requirement protein A2, OMI, HtrA2, PRSS25, Serine protease 25, Serine proteinase OMI, Serine protease HTRA2, mitochondrial, HTRA2, Omi stress-regulated endoprotease

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.