



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Phospho-KCNB1 (pS563) Monoclonal Antibody, L83/11R

OVERVIEW

PRODUCT REFERENCE	PTX19239
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	h-DRK1, Delayed rectifier potassium channel 1, Potassium voltage-gated channel subfamily B member 1, KCNB1, Voltage-gated potassium channel subunit Kv2.1, DRK1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.