



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Phospho-PLC- γ (pY783) Monoclonal Antibody, SAA0303

OVERVIEW

PRODUCT REFERENCE	PTX19241
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Phosphoinositide phospholipase C-gamma-1, PLC1, PLC-II, Phospholipase C-II, Phospholipase C-gamma-1, PLC-gamma-1, PLC-148, PLCG1, 1-phosphatidylinositol 4, 5-bisphosphate phosphodiesterase gamma-1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human, Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.