



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Phospho-CASC5 (pT943) Monoclonal Antibody, pMELT-KNL1

OVERVIEW

PRODUCT REFERENCE	PTX19242
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Blinkin, CT29, Kinetochores-null protein 1, Cancer/testis antigen 29, Bub-linking kinetochores protein, Protein CASC5, Protein D40/AF15q14, ALL1-fused gene from chromosome 15q14 protein, Kinetochores scaffold 1, CASC5, KNL1, Cancer susceptibility candidate g

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.