



## Anti-Human TSLPR/CRLF2 Monoclonal Antibody, 13H5

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19243
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	Cytokine receptor-like 2, CRLF2, CRL2, Thymic stromal lymphopoietin protein receptor, Cytokine receptor-like factor 2, TSLPR, ILXR, IL-XR, TSLP receptor

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

