



Anti-Human TIMP1 Monoclonal Antibody, 2PC10

OVERVIEW

PRODUCT REFERENCE	PTX19246
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Erythroid-potentiating activity, TIMP, EPA, TIMP1, Metalloproteinase inhibitor 1, TIMP-1, Fibroblast collagenase inhibitor, Collagenase inhibitor, CLGI, Tissue inhibitor of metalloproteinases 1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

