



Anti-Human TIMP2 Monoclonal Antibody, SAA0317

OVERVIEW

PRODUCT REFERENCE	PTX19247
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	TIMP2, TIMP-2, CSC-21K, Tissue inhibitor of metalloproteinases 2, Metalloproteinase inhibitor 2

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

